

PEX11B-peroxisomal biogenesis factor 11 beta Gene

Reference: PTXBC011963

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PEX11-BETA; PEX14B; peroxisomal membrane protein 11B; peroxin-11B; peroxisomal biogenesis factor 11B; protein PEX11 homolog beta; peroxisomal biogenesis factor 11 beta		

For research use only.

GENE DESCRIPTION	peroxisomal biogenesis factor 11 beta
GENE ID	8799
ORIGIN SPECIES	Human
ACCESSIONS	BC011963
VECTOR	pDONR223
NCBI SYMBOL	PEX11B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgctgggtccgcttcagtgctcagagccaagccgggagcggctgtgtagggccgcccagatgcttctctcttggccatgcgctgcagaggcagtgagccagctgagcttggaaagaagcttctacgc
ctgggttaactcagcagatgccttgagtcagccaaagagctgttcacatcagatgttgcctgagattcgcacacactgttagtcacacacacagcctgtgactctgtgtgacaatgtcctgtggcctggaaagtctggactggctccccgtgtgatcaggagaaagtggccagcg
ttcattcaggtactattgtttccctcatcatgaatttgagccgtgatcctatgagattgcctactgatggagcaagagctcttctgtgagccggcgactgaaaggtctgtggaggaggatcccaaggaggaagtgaactggggacttggggaccaggagactcaggaggaggtctg
ccccaaactggctctgaaactcggctgcaagtctgtcctggctcagtccttagaggtcatccccactctgttagacgtggtcagaaatgctgtgatcttcatctctggacaaactaggcctctggcctgtggcctgggattgtgggcttggcctctgtctccatctgt
ctattctaccctaatactatccctggctacgactcaagccctga