

FAM100A-family with sequence similarity 100, member A Gene

Reference: PTXBC012018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM100A; FAM100A; PP11303; UBA-like domain-containing protein 1; family with sequence similarity 100, member A; UBA like domain containing 1		

For research use only.

GENE DESCRIPTION	family with sequence similarity 100, member A
GENE ID	124402
ORIGIN SPECIES	Human
ACCESSIONS	BC012018
VECTOR	pDONR223
NCBI SYMBOL	FAM100A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgtgaacatggacgagctcaagcaccaggtcatgataaccagttcgtgtgacggcggtgctgcggcgccgaccaggcgaagcaactgctgcaggcggcccaactggcagttcgagacagccctcagcgcttttccaggagaccaacatccctacagccaccatcaccac
cagatgatgtgacccccccaataccctgtacaccccccaacttccctgacgtctcaccatgttctccgtctcaaggcctccgagagcttccacagcgggtggcagcggcagcccgatggccgcgacagccacgtcaccctccacacttccccatgccgcccagcagctc
tgcggcctccagctggcccacggcggtcgtgccccggggggccacagcaccaccagccagccgcccctgtggactccaaccccccttctcggcttcagactggccacccctggcccccaacaggccacctcagaacccagggcccccctgcatggagcagagag
ataa