

NDUFS1-NADH dehydrogenase (ubiquinone) Fe-S protein 1, 75kDa (NADH-coenzyme Q reductase) Gene

Reference: PTXBC012068

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CI-75Kd; CI-75k; PRO1304; NADH-ubiquinone oxidoreductase 75 kDa subunit, mitochondrial; NADH dehydrogenase (ubiquinone) Fe-S protein 1, 75kDa (NADH-coenzyme Q reductase); complex I 75kDa subunit; complex I, mitochondrial respiratory chain, 75-kD subunit;		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) Fe-S protein 1, 75kDa (NADH-coenzyme Q reductase)
GENE ID	4719
ORIGIN SPECIES	Human
ACCESSIONS	BC012068
VECTOR	pDONR223
NCBI SYMBOL	NDUFS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggttttaggcagttctgcactccaaagaatgatggagcagcaattctgcagctgtttctagcattgcacaaagattcgatgactagtgtgttactggtgattggaaagtatgaatatcttcataaggattgcaagtaagtagctgtttggaccttggtctataagcctgggggtga
agcaattcggagaacccctccaaaggtgctgtttctcctgggagcagatggaggtgtatcacacgacagagattgccaaagattgttcattattatcaaggacatcatggtgatgttgggctcccatagctgatgttattctccaggagctgttacacagagaagctctgctacatatgc
aacactgagggtagagctcagcagactaaggtagcagtgacacctctgcttggcaagagaagactggaaaattataagagcactctctgagattgctggaatgactctccatgatactctggatcaagtaaggaaacagattggaagaagctctctcctaattctgtcgatatgatgatt
gaaggggctaattacttcagcaagcaaatgagctctcaaaagctagtgaaccagcagcttctgtgaccacctgtttccacctgactataaaagacttctacatgacagattcaattagcagagcctcacagacaatggccaaatgtgtcaaaagctgtcacagagggtgccaggca
gtagaggaacctcatatgctga