

RABGAP1L-RAB GTPase activating protein 1-like Gene

Reference: PTXBC012094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HHL; TBC1D18; rab GTPase-activating protein 1-like; TBC1 domain family, member 18; expressed in hematopoietic cells, heart, liver (HLL); RAB GTPase activating protein 1 like		

For research use only.

GENE DESCRIPTION	RAB GTPase activating protein 1-like
GENE ID	9910
ORIGIN SPECIES	Human
ACCESSIONS	BC012094
VECTOR	pDONR223
NCBI SYMBOL	RABGAP1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggaagaatctccattatgtagcctatgacgcccatgtttcagccagctgcacgatgaagacttcctcactagtctgtggccatcagcaagccaggtctatggtaccaaccaagaagctgaagaaatagagaaagaatatcagacaatgcgagagagtcagctgcaacaggaagacccaatggatagatacaagtttgtatattgttag