



Reference: PTXBC012100

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	integral membrane protein GPR137B; transmembrane 7 superfamily member 1 (upregulated in kidney); transmembrane 7 superfamily member 1 protein; G protein-coupled receptor 137B		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 137B
GENE ID	7107
ORIGIN SPECIES	Human
ACCESSIONS	BC012100
VECTOR	pDONR223
NCBI SYMBOL	GPR137B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]