



Reference: PTXBC012102

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BNAH2; receptor-type tyrosine-protein phosphatase F; LCA-homolog; leukocyte antigen-related (LAR) PTP receptor; leukocyte antigen-related tyrosine phosphatase; leukocyte common antigen related; protein tyrosine phosphatase, receptor type, F polypeptide; r		

GENE DESCRIPTION	protein tyrosine phosphatase, receptor type, F
GENE ID	5792
ORIGIN SPECIES	Human
ACCESSIONS	BC012102
VECTOR	pDONR223
NCBI SYMBOL	PTPRF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1