

LOC440700-carbonic anhydrase XIV (CA14) pseudogene Gene

Reference: PTXBC012108

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	carbonic anhydrase XIV (CA14) pseudogene
GENE ID	440700
ORIGIN SPECIES	Human
ACCESSIONS	BC012108
VECTOR	pDONR223
NCBI SYMBOL	LOC440700
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgggtggctcacactgtatcccagcacttgggagggcggaggtgggtggattgcctgagctcaggagttccagaccagcctggacaacatgggtgacttccagaaaatgtagctaccagctacacctgcactggggtcaggatga