



Reference: PTXBC012121

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APNH; LIKNS; NHE-1; NHE1; PPP1R143; sodium/hydrogen exchanger 1; Na(+)/H(+) exchanger 1; Na-Li countertransporter; protein phosphatase 1, regulatory subunit 143; solute carrier family 9 (sodium/hydrogen exchanger), isoform 1 (antiporter, Na+/H+, amiloride		

For research use only.

GENE DESCRIPTION	solute carrier family 9 (sodium/hydrogen exchanger), member 1
GENE ID	6548
ORIGIN SPECIES	Human
ACCESSIONS	BC012121
VECTOR	pDONR223
NCBI SYMBOL	SLC9A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr