

TMEM106A-transmembrane protein 106A Gene

Reference: PTXBC012139

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 106A		
For research use only.			

GENE DESCRIPTION	transmembrane protein 106A
GENE ID	113277
ORIGIN SPECIES	Human
ACCESSIONS	BC012139
VECTOR	pDONR223
NCBI SYMBOL	TMEM106A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtaagacgttttccagctgggctcttggcgggaggtatgagaacaagtcaatctgtctcctcaaacagccattggcagcaaggctgtcaactaccagcacccgttagcagcaagctttttgttctgtgtgcctgtgaaggaaactgctgatgccagcttctgactgtccacctg
ccaggcgagtggaagattcccaagagctggagaagcagttgggtgctctatccctatggggaccagaggtgaagcccaagcacacgaagctctttgttctgtgcccgtgtcctatctgctgtgacctctctctatcgtcttttctgtttcccggtccgtcattgtgcagctg
caggcctcaactcctccagctggcctttgatgaggctgatactacacataacgaatacttaaacatctccaatggcaactactaccctattgtgtgacacagctgacctcgagggttctgcacctgtccctcgtgtgtgggaggtttccaacaaccttctctacacattggccctt
ggccagtgaacagatgttttacgcagtagctaccaagatagggatgaaacacatacaaaatctgtacctgggtggaaatcaaatgcccattgtctttgcacatccaggccacctgacctgttcatctgagccattcagagcagctggtcttcagagctatgaatatgtggactgc
cgaggaaacgcactgtgtcccccagctgacctcaccacattga