



Reference: PTXBC012142

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6C2; VMA5; V-type proton ATPase subunit C 2; ATPase, H ⁺ transporting, lysosomal 42kDa, V1 subunit C2; V-ATPase C2 subunit; vacuolar proton pump subunit C 2; ATPase H ⁺ transporting V1 subunit C2		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal 42kDa, V1 subunit C2
GENE ID	245973
ORIGIN SPECIES	Human
ACCESSIONS	BC012142
VECTOR	pDONR223
NCBI SYMBOL	ATP6V1C2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr