

SNX27-sorting nexin family member 27 Gene

Reference: PTXBC012184

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRT1; MY014; sorting nexin-27; methamphetamine-responsive transcript 1; sorting nexin family member 27		

For research use only.

GENE DESCRIPTION	sorting nexin family member 27
GENE ID	81609
ORIGIN SPECIES	Human
ACCESSIONS	BC012184
VECTOR	pDONR223
NCBI SYMBOL	SNX27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagtacacagtgaattactttgccttatttgaagtgcagtcactcctttgtacgtaaattggcacctaatgagtttctcacaaactctacattcagaattatatacagctgtgccaggccactgcttgaccattcgaagtggcttttacaacagaagaagaattctttaaagaca
atgacctgtctgttacctacttctttcatcaggcagtcgatggtgaagaagggttacatcaaaagcagaagaaaaagtcctatcaattacagaagctatacgaacaaagaaaatggctacgtacccaacatgctaaggactgtgagggtcacaatgaatcatcttccccactgtgcctgtga
ctccaggaggaaagggcagctttatcacagccatcagcatcacgcactttaaactgcatgcctgcactgaagaaggacagctggagaaccaggtaattgcaattgaatgggatgagatgcagcgtgggacacagatgaagaaggatggccttctgttcgaatgacagaggagaga
agaagccccgatgggttaaaatcttcacgccatatttcaattacatgcatgagtgcttcgagaggggtgtctgcgagctcaagtggaagaaagaggaaatattag