

C20orf27-chromosome 20 open reading frame 27 Gene

Reference: PTXBC012196

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0687 protein C20orf27; chromosome 20 open reading frame 27		
For research use only.			

GENE DESCRIPTION	chromosome 20 open reading frame 27
GENE ID	54976
ORIGIN SPECIES	Human
ACCESSIONS	BC012196
VECTOR	pDONR223
NCBI SYMBOL	C20orf27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgagccaagaaggcaagaagccagagtcggagatccgcttgcggcagggccacgatgcagaaggatcccacagccagtccttgcagagaagctgcagtcgggtcatggtcaccagggagagtgacagcagcttctggtaaggtggcttctgaagat
cctgcacaggtatgagattactctactctgccccagtgacagggctgagcaaggatgtccgcgagggcactgtccccagcctgcacctaagctctcagcgtggtgcccgtccctgaaggttatagtgtaagtgtagtactggcgacaaagagggcgtcctcaagaggagat
actgctagcctgcgaaggtggcactggcacctgtgtgcgctgacgggtgacggcccgctcatggaccggcaccacggcagcccatgctgctggatggtgtcaagtggtggcgccgagctggaatagactcagagcagcagcactggcacggcttgactga