

C11orf54-chromosome 11 open reading frame 54 Gene

Reference: PTXBC012298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ester hydrolase C11orf54; PTD012; PTOD012; chromosome 11 open reading frame 54		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 54
GENE ID	28970
ORIGIN SPECIES	Human
ACCESSIONS	BC012298
VECTOR	pDONR223
NCBI SYMBOL	C11orf54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtgctgagttttttcatgtaccaagcttgaagagcttgcggagttatgcagaaggggttaaaagataacttgcgtgatgccaggtctctgtagttgattgccctgatttgactaaggaaccctttaccttctgtaaaaggcatctgtgggaaaactagaattgcagaagttggag
gtgtgcttacttattgcctctgttaaaccaaaaaaatttatgatctgaataaaatgcagaagaaatcaagctgcctggagcctttattcttggagcaggagcagggtccatttcagactctcgggttcaattctgagtttatgccagttatcagacagaagtgacacaaagcctctgtaaatg
gaagttactttgccatgtgaaccctgcagatggagggtgcctactggagaatacagtgagaatgtcatgatttcagtgctcattactggcctaactttttgccagtgaaagccaactggcaagcctgcagaattttctcctgccccctgaactctgatgaagaagtgaataaatgtgtgcat
ttttatgaaatgaagctcctttgtgtgtctaccagttttgtctccagagaccaggggttgatttgcgactggagcacactcatttttttagtctcatggagaaggtggacactaccattatgacactactccagatatagtggaatattcttggaatactttacctgcagagtttctctatgcattga
tcaacaaaaagagacgattcaattggcgagattaa