

FABP7-fatty acid binding protein 7, brain Gene

Reference: PTXBC012299

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B-FABP; BLBP; FABPB; MRG; fatty acid-binding protein, brain; brain lipid-binding protein; brain-type fatty acid-binding protein; mammary-derived growth inhibitor-related; fatty acid binding protein 7		

For research use only.

GENE DESCRIPTION	fatty acid binding protein 7, brain
GENE ID	2173
ORIGIN SPECIES	Human
ACCESSIONS	BC012299
VECTOR	pDONR223
NCBI SYMBOL	FABP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtggaggcttctgtgctactggaagctgaccaacagtcagaacttgatgagtacatgaaggctctagggctggccttgcactagcaggtgggaaatgtgaccaaacaacggtaattatcagtcagaaggagacaaagtggatcaggactctcagcacattcaagaacac
ggagattagttccagctgggagaagatttgatgaaccactgcagatgatagaactgtaagtctgtgttagcctggatggagacaaactgttcacatacagaatggatggcaagaacaattttgaagagaataaggatggcaaatggttatgaccttacttttggatgt
ggttgctgttcgacctatgagaaggcataa