

## PDZK1IP1-PDZK1 interacting protein 1 Gene

Reference: PTXBC012303

|                |                                                                                                                                                                                                                               |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                           | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                          | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | DD96; MAP17; SPAP; PDZK1-interacting protein 1; 17 kDa membrane-associated protein; epithelial protein up-regulated in carcinoma, membrane associated protein 17; membrane-associated protein 17; PDZK1 interacting protein 1 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | PDZK1 interacting protein 1                                                                       |
| <b>GENE ID</b>          | 10158                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC012303                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PDZK1IP1                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtcgccctcagcctctctctggtgctcatggcagtgccacctgccagctgcagcaaggcctgggaaccttcagccctggatgcaggcccttatcgcggtggcgtgttctgtgtctcgttgcaatcgcccttcagtcgaacactctctgtgccaggaggagccggagcc  
tgcacacatgatcctgaccgtcggaaacaaggcagatggagtcctgtgtgggaacagatggaaggctactcttcgatggcggccagtttcagggtccagtgagcatgagaatgcctatgagaatgtgcccgaggagggaaggcaaggtccgagcaccgccgatgtaa