

LINCR-likely ortholog of mouse lung-inducible Neutralized-related C3HC4 RING domain protein Gene

Reference: PTXBC012317

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LINCR; E3 ubiquitin-protein ligase NEURL3; lung-inducible neuralized-related C3CH4 RING domain protein; neuralized homolog 3 pseudogene; neuralized-like protein 3; neuralized E3 ubiquitin protein ligase 3		

For research use only.

GENE DESCRIPTION	likely ortholog of mouse lung-inducible Neutralized-related C3HC4 RING domain protein
GENE ID	93082
ORIGIN SPECIES	Human
ACCESSIONS	BC012317
VECTOR	pDONR223
NCBI SYMBOL	LINCR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgagccaaggcgcacaggtcgtctggacacgcgtggctgcatcgcgcacaggcaccacgttccacgagcgcctggttcagccagcggcggcgcctggcgagcgttggcctgcgagtgctcgggagagagcggctggtcggcggcctccgctg
ggcttcacggcctggacccccgctgctgctccgctgccagcctgccgccctcctgctgccccgacctggaggagcagagcccacgtggcgccctgctgctgagggctgcgcctcactgggacttggctccgctctgggtgaccgccggcctctcccaaggtc
aacgccggctgccgctcctgctgctgagggcgtgccctcggcgccccgctctgggcctga