



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

SLC22A5-solute carrier family 22 (organic cation/carnitine transporter), member 5 Gene

Reference: PTXBC012325

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDSP; solute carrier family 22 member 5; high-affinity sodium dependent carnitine cotransporter; organic cation/carnitine transporter 2		

For research use only.

GENE DESCRIPTION	solute carrier family 22 (organic cation/carnitine transporter), member 5
GENE ID	6584
ORIGIN SPECIES	Human
ACCESSIONS	BC012325
VECTOR	pDONR223
NCBI SYMBOL	SLC22A5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr