

DNAJA3-DnaJ (Hsp40) homolog, subfamily A, member 3 Gene

Reference: PTXBC012343

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCA57; TID1; hTID-1; dnaJ homolog subfamily A member 3, mitochondrial; DnaJ (Hsp40) homolog, subfamily A, member 3; dnaJ protein Tid-1; hepatocellular carcinoma-associated antigen 57; tumorous imaginal discs protein Tid56 homolog; DnaJ heat shock protein		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily A, member 3
GENE ID	9093
ORIGIN SPECIES	Human
ACCESSIONS	BC012343
VECTOR	pDONR223
NCBI SYMBOL	DNAJA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttgacattcaatcaagctgcaaaaggggtcaacaaggagttcaccgtgaacatcatggacacgtgtgagcgtgcaacgggcaagggaacgagccggcaccagggtgcagcattgccactactgtggcggctccggcatggaaacatcaacacaggcccttttgatgc
gttccacgtgtaggagatgtgtggcggcgtccatcatcatatcgcctgtgtgtgtcaggggagcaggacaagcagaagaaagcagtgatgatccctgtgcctgcaggagtcgaggatggccagaccgtgaggatgcctgtgggaaaaaggaaatttcattacgttca
gggtgcagaaaagccctgtgtccggaggacggcgagacatccactccgacctttatttctatagctcaggtctctctgggggaacagccagagcccggtgtacgagacgaacgtgacgatccccctgggactcagacagaccagaagattcggatgggtggaaa
ggcatccccggattaacagctacggctacggagaccactatccatcaagatagcttcaagaggctaacgagccggcagcagagcctgatcctgagctacccgaggacgagacagatgtggaggggacgtgaacggcgctcacctcaccagctctgaaaaagatc
cactggaactag