



Reference: PTXBC012355

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	creatine transporter SLC6A8 variant D; CCDS1; CRT; CRTR; CT1; CTR5; sodium- and chloride-dependent creatine transporter 1; creatine transporter 1; solute carrier family 6 (neurotransmitter transporter), member 8; solute carrier family 6 (neurotransmitter		

For research use only.

GENE DESCRIPTION	solute carrier family 6 (neurotransmitter transporter, creatine), member 8
GENE ID	6535
ORIGIN SPECIES	Human
ACCESSIONS	BC012355
VECTOR	pDONR223
NCBI SYMBOL	SLC6A8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]