

KLC4-kinesin light chain 4 Gene

Reference: PTXBC012357

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KNLSL8; bA387M24.3; kinesin light chain 4; kinesin-like protein 8		
For research use only.			

GENE DESCRIPTION	kinesin light chain 4
GENE ID	89953
ORIGIN SPECIES	Human
ACCESSIONS	BC012357
VECTOR	pDONR223
NCBI SYMBOL	KLC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaggcctgtgttggggcagcgggatgagcctgcagggccacggctcagccaagaggagatctggggagcacagcgtgtcagccaaggctagaggccctacgcagtgaaacaccaggccgtgtgcaaaagcctgtccagaccattgagtgtcagcaggaggcc
atgaggaaaggctgtgcatgagaagcccgagcttccgcttctatggaaaacattgagctcgggctgagtgaggccaggtgatgctggtctagccagccactgagcacagtgaggctggagaaacagaagctgcgggctcaggtgcgggctatgccaggagaaccag
tggtctcgggatgagctggctggcaccagcagcggctacagcgcagtgaaacaggctgtggctcagctggaggaggaaaagaagcacctggagttcctggggcagctgcggcagtgatgaggatggacatacctcggaggagaaagagcgatgccaccaaggattccctgg
atgaccttttctaagaggaggaagaggacccagcaatggctgtccctgtgctcaaggtgtctacagcagctcagcagggtggatagatgccagcaaggtgtgcgacgtgtgcacaacctgtgtatccagtcagcagccaaggtgtgctatgagtggtggcgtgccactctgtaagc
aggcactagaggacctggagcgcacatcagccgtggccacctgatgtcgccacatgctcaacatccttcttgggtatcgtgaccagaataagtaaggaagctgccacctgtgaatgatgcccttagcatccggagagacacttgggacctgacctcctgtgtcagtatt
ccttgcctccccaccccgcccgacccccaccattgctgttttggcctcttag