

DLX3-distal-less homeobox 3 Gene

Reference: PTXBC012361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AI4; TDO; homeobox protein DLX-3; distal-less homeobox 3		
For research use only.			

GENE DESCRIPTION	distal-less homeobox 3
GENE ID	1747
ORIGIN SPECIES	Human
ACCESSIONS	BC012361
VECTOR	pDONR223
NCBI SYMBOL	DLX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtggtctcttcgatgcgaagctcagcagcatctcaccgacatctccagctcccttagctgccatgcgggctccaaggactcgcctaccctgccgagcttctgtcactgacctgggctactacagcgtccccagcagcattactactcgggcagccctatggccagcgggtgaa
ccccctacacctaccaccaccaatcaatcgaatgggcttcagggcagcgggcgttactcgcctcgaagtcggaatatactacggagcctcctaccggcaatacggggcgtatcgggagcagccgctgccagcccaggaccagtgctgggtaaggaggagccggaagcagagtgctc
gcatggtgaatgggaagcccaagaaggtccgaaagccggtacaatctactccagctaccagctgcccgcctgcagcgccttcagaaagccagctacgtggcgtgcccagcgcgcgcagctgcccgcgcagctgggcctcagcagacacaggtgaaatctggttcca
gaaccgcgttccaagttcaagaactctacaagaacggggaggtgcccgtggagcagtcctcaataacagtgattccatggcctgcgaactcaccacatcacccgcctctgggacacctctccactccactcggccctgcccgcagtcagctgccccgcctccataca
gtgcctccccagctacctggagacccccaccaactcgtggtatcacgcagaaactgagtggaacccacttacagcagcagccctcagccagccaccctgcacatgcctctccgggcccccccaacctgggctgtgtactga