

LBH-limb bud and heart development homolog (mouse) Gene

Reference: PTXBC012373

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein LBH; hLBH; limb bud and heart development homolog; limb bud and heart development		

For research use only.

GENE DESCRIPTION	limb bud and heart development homolog (mouse)
GENE ID	81606
ORIGIN SPECIES	Human
ACCESSIONS	BC012373
VECTOR	pDONR223
NCBI SYMBOL	LBH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggggcaggaggacgaggggtgcgtgtggacacagcagtcgcgcgaattcgtttctgggaagccaatgtcgcggcacccttgccttcctctgtgtgtgctgtgtgacacacatcaatggcaatacttctccaactcctcgcagaagtgggagagggcggcagcctg
caccgagaggggctttctctctctgtcctccgccttcgtttgtgtgctgcagagagtggttcacatactctcattccctcgcctcccttgtggacgggggtcttgcctttcaattctgtgtttgtgtcttccctatctgctaccctgaatcacctgtcctgtctgtgtgatgggaa
catgcttataa