

TSPAN6-tetraspanin 6 Gene

Reference: PTXBC012389

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T245; TSPAN-6; tetraspanin-6; A15 homolog; tetraspan TM4SF; tetraspanin TM4-D; transmembrane 4 superfamily member 6; tetraspanin 6		

For research use only.

GENE DESCRIPTION	tetraspanin 6
GENE ID	7105
ORIGIN SPECIES	Human
ACCESSIONS	BC012389
VECTOR	pDONR223
NCBI SYMBOL	TSPAN6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggcgtcccgtctcggagactgcagactaaccagtcactactgtttcaagagcgttctgctaactctacactttatcttgatcactggcgttctcttcagttggcatttggcaagggtgagcctggagaattactttctcttttaaatgagaagccaccaatgtcccctcgtgctc
attgctactggtaccgtcattatcttttggcacccttgggtgtttgctacctccgagcctctgcatgcatgctaaaactgtatgcaatgtttctgactctcgttttttggctgaaactggtcgtcgtccatcgtaggattgttttcagacatgagattaagaacagctttaagaataatgaaaaagct
ttgaagcagataaactacaggagattatagaagccatgcagtagacaagatccaaaatcgttgcattgttgggtgtcaccgattatagagattggacagataactaattactcagaaaaaggatttcttaagattgctgtaaaacttgaagattgtactccacagagagatgcagacaa
gtaaacatgaaggttgtttataaaggatgaccattatagagtcagaaatggagtcgttcaggaattcttggagtgcttccaactgattggaatcttctcgtcctctcgtgcataacaataaccagatgagatagtgtaa
```