



Reference: PTXBC012418

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACAT; CE-1; CEH; CES2; HMSE; HMSE1; PCE-1; REH; SES1; TGH; hCE-1; liver carboxylesterase 1; acyl coenzyme A:cholesterol acyltransferase; brain carboxylesterase hBr1; carboxylesterase 1 (monocyte/macrophage serine esterase 1); carboxylesterase 2 (liver); c		

For research use only.

GENE DESCRIPTION	carboxylesterase 1 (monocyte/macrophage serine esterase 1)
GENE ID	1066
ORIGIN SPECIES	Human
ACCESSIONS	BC012418
VECTOR	pDONR223
NCBI SYMBOL	CES1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr