

RND3-Rho family GTPase 3 Gene

Reference: PTXBC012513

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARHE; Rho8; RhoE; memB; rho-related GTP-binding protein RhoE; protein MemB; ras homolog gene family, member E; rho-related GTP-binding protein Rho8; small GTP binding protein Rho8; Rho family GTPase 3		

For research use only.

GENE DESCRIPTION	Rho family GTPase 3
GENE ID	390
ORIGIN SPECIES	Human
ACCESSIONS	BC012513
VECTOR	pDONR223
NCBI SYMBOL	RND3
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggagagaagccagccagaaatccagcaaatctatcatggatcctaatacagaactgaaatgcaagatagttggtgggagacagtcagtggaataactgctgctccatgtctcccaagactgctccccgagaattacgttctcagtggttgaattacag
gcccagtttgaatcgacacacaagaatagagtgagcctgfgggacactcgggttctctactatgacaatgcccggcctctcttacctgattcggatgctgtgctgattgctttgacatcagtagaccagagaccctggacagtgctcctcaaaagtggaaagtgaaatccagga
atfttgccaaataccaaaatgctcttggcggctgcaagtctgctcggacagatgtagtacattagtagagctctccaatcacaggcagacgccagtgctctatgaccaggggcaaatatggccaacagattggagcagctactatcgaatgctcagcttacagtcgaaaata
gcgtcagagacattttcagcttggccacctggcatgftaataagacaataaaaacgtaagcgggaacaaatcacagagagccacaagcggattcaccatgcctagcagaccagaactctcggcagttgctacggacttacgaaagacaagcgaagagctcactgtgatgtg
a