

ZNF673-zinc finger family member 673 Gene

Reference: PTXBC012569

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF673; KRAB domain-containing protein 4; KRAB box domain-containing protein 4; zinc finger family member 673; zinc finger protein 673; KRAB box domain containing 4		

For research use only.

GENE DESCRIPTION	zinc finger family member 673
GENE ID	55634
ORIGIN SPECIES	Human
ACCESSIONS	BC012569
VECTOR	pDONR223
NCBI SYMBOL	ZNF673
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccatgtccaggaaatcattgacctcaaggacgtgtttgtggacttcacctggaggagtgccagcaactggactctgccagaagaacctctacagggatgtcatgcttgagaactacagccacctgggtccgtgggtatctagttgcgaagccagatgtgatctcaggttggga
ccaggtgaagagtcctggatggcagatgggggaccccggtacggacctgtgcaggtgaggacagggccagatgtatccattttgcccctatgtatttgaagtgcgtgtattag