

ZC3H7A-zinc finger CCCH-type containing 7A Gene

Reference: PTXBC012575

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC055; ZC3H7; ZC3HDC7; zinc finger CCCH domain-containing protein 7A; zinc finger CCCH-type domain containing 7; zinc-finger protein AY163807; zinc finger CCCH-type containing 7A		

For research use only.

GENE DESCRIPTION	zinc finger CCCH-type containing 7A
GENE ID	29066
ORIGIN SPECIES	Human
ACCESSIONS	BC012575
VECTOR	pDONR223
NCBI SYMBOL	ZC3H7A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggagaatgggatacaagataggagcaattttacgaactatggctcaagagtgcaaaaaatgaaaaagtgaagacatagccagtcagtcagcaacaaggaaaatggaaaacaaatccatgccacagattatgctgaagttacagtggactttcactgctggatgtgtgggaaaaa
ctgcaacagtgagaagcagtgccaggccacatctctccgagaagcacaagagaagggtttccacaccgaggagcagcagctactgctggcagcaccgcttcccaacaggctatttcagtatttggatagggtatgaatggcactgccagaggaaacagctgaaattgacat
ggaaatgccgaacttcatgaatgggaagaagaagagatgccctaaagatgaagctcaacaagcagcaaaagatcacttaattggcccaatgataatgacttggaaaatagttttttgttaaagatttaaactaa