



Reference: PTXBC012606

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOV34; P40; Rpn8; S12; 26S proteasome non-ATPase regulatory subunit 7; 26S proteasome regulatory subunit S12; 26S proteasome regulatory subunit rpn8; Moloney leukemia virus-34 proviral integration; Mov34 homolog; proteasome (prosome, macropain) 26S subuni		

GENE DESCRIPTION	proteasome (prosome, macropain) 26S subunit, non-ATPase, 7
GENE ID	5713
ORIGIN SPECIES	Human
ACCESSIONS	BC012606
VECTOR	pDONR223
NCBI SYMBOL	PSMD7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1