

GUCY1A3-guanylate cyclase 1, soluble, alpha 3 Gene

Reference: PTXBC012627

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GC-SA3; GUC1A3; GUCA3; GUCA3A; GUCY1A1; MYMY6; guanylate cyclase soluble subunit alpha-3; GC-S-alpha-1; GCS-alpha-3; guanylate cyclase 1, soluble, alpha 3; soluble guanylate cyclase large subunit; guanylate cyclase 1 soluble subunit alpha		

For research use only.

GENE DESCRIPTION	guanylate cyclase 1, soluble, alpha 3
GENE ID	2982
ORIGIN SPECIES	Human
ACCESSIONS	BC012627
VECTOR	pDONR223
NCBI SYMBOL	GUCY1A3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttctgcacgaagctcaaggatctcaagatcacaggagagtgctcttctcttactggcaccaggtcaagttcttaacgagcttccagaggcagcaggagctcgaagcaaccgtgcccatctgtcaagacattcctgagaagaacatacaagaagtcttctca
aagaaaaaccagtcggagccgagctatcttcacacttggcagagagatatttgcacactgatttcccagagtttgaacggctgaatgttgcacttcagagaacattggcagaacacataaaagaaagcaggaaatcttggaaagagaagactttgaaaaacaattgcagagcaagc
agttgcagcaggtggagaccattggcgatgcctattgttagctggggattacacaaagagagtga