



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CHMP1B-chromatin modifying protein 1B Gene

Reference: PTXBC012733

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf2; C18-ORF2; C18orf2; CHMP1.5; Vps46-2; Vps46B; hVps46-2; charged multivesicular body protein 1b; chromatin modifying protein 1B; chromatin-modifying protein 1b; vacuolar protein sorting 46-2; vacuolar protein sorting-associated protein 46-2		

For research use only.

GENE DESCRIPTION	chromatin modifying protein 1B
GENE ID	57132
ORIGIN SPECIES	Human
ACCESSIONS	BC012733
VECTOR	pDONR223
NCBI SYMBOL	CHMP1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctaacatggagaacacacctgttcaacctgaagttcgcggccaaagaactgagtaggagtgccaaaaaattcgataaaggaggaaaaggccgaaaaaggccaaaattaaaaaggccattcagaaggggcaacatggaaagtgcggaggatacacgccaaaaatcccatccgcagaag
aacacggcgggtgaattttctgagaatgagtgccgacgtcgtatgcagtggttcgccagggtccagcagccggctgacgatggggcaagtgtagcaaatgctgatggctgtgttgggttaatgctgatggatgcacattgaagacacatgaatctggagagaatttctccttgatggacaatttcagc
accagttgagctactctggacgtccagacgaccaaattggaagacacgatgagcagcagcagcagctccacactcccgacaacaaatggatattgctctcgagaattgagctgctgcgcagatgaggcgggctcgacctcaacatggagctgcgcaggccagacagcgctcctgtggccac
gacgtctggatcgcggacacgaattgaactctcagaagctgccccttcgcgtatctgaa

More info: Please contact us for any question: contact@proteogenix.fr