



Reference: PTXBC012759

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZFP64 zinc finger protein; ZNF338; zinc finger protein 64; zinc finger protein 338; zinc finger protein 64 homolog		

For research use only.

GENE DESCRIPTION	zinc finger protein 64 homolog (mouse)
GENE ID	55734
ORIGIN SPECIES	Human
ACCESSIONS	BC012759
VECTOR	pDONR223
NCBI SYMBOL	ZFP64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]