



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

XRCC6BP1-XRCC6 binding protein 1 Gene

Reference: PTXBC012776

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	XRCC6BP1; mitochondrial inner membrane protease ATP23 homolog; Ku70-binding protein 3; XRCC6-binding protein 1; ATP23 metallopeptidase and ATP synthase assembly factor homolog		

For research use only.

GENE DESCRIPTION	XRCC6 binding protein 1
GENE ID	91419
ORIGIN SPECIES	Human
ACCESSIONS	BC012776
VECTOR	pDONR223
NCBI SYMBOL	XRCC6BP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggcggagctccgcagagcgcggcggggccggccgcagggggagcagctgcagcagcaacagctctctccaggtctcccgagcgtctggccacggggaaatcccgacaaagggtctctcagctcttcacagcaacagaagatgccagcttaggtctctgaagac
gtcggagacaaatccatcatgtccaaactctctctgctgagctctgaacacactcagggtgtgctgttaacaagaatagacacactcttctcggaaagcgttaatggaaatgtcagctggaggttttgatctctcaactcagatgtttgttggccagaataaataatcataatgccgcccataatgaacagagt
gtgacacacagcagctattaccctcatttgatcattatctgcgcggagctgcagctgttgcacacagacatttggcgctgcagctgctgcagctgctcaacttctggagacagctctcactctcaatgaataatcaggtctacatttggtaataaacacacacccacccagcactgtgtgctgcagaca
gaccactctctctatctctgtctgtaggaatcagcacaagaagaatgtctaaagaagcgtgtgtgaagtgtttgaaactctgtctaaagaccatgaaccttctgttggaaagctccacataaacagatctcaagaatgtctcacaagagatctgaaaacctgtacgtgtaattctcaaatataga

More info: Please contact us for any question: contact@proteogenix.fr