

FOXR2-forkhead box R2 Gene

Reference: PTXBC012934

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	forkhead box protein R2; forkhead box protein N6; forkhead box R2		
For research use only.			

GENE DESCRIPTION	forkhead box R2
GENE ID	139628
ORIGIN SPECIES	Human
ACCESSIONS	BC012934
VECTOR	pDONR223
NCBI SYMBOL	FOXR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactaaactaaagactgtgaattttggtatagctccatggccaggtcccagggtgctgactgggacatgaggatgagttatttctgctgtaccacagaccagtgtcttttagctgagcaaatccttgccaaatacagagtcggagtaataagccccagaaatgcctcaga
agaggagaccagtcctgatggagatggctccctgtgaaccaatctgtggtatgtgggtgaccccaatctgtgcccccttggcagccaggaggccccaaagcccagtggaagaggatctgacaaatttctcttccctcagcccccaaaaagacgaagggtcctaact
gtcagaggacaaagtgttagagctctgccccttccagtgagcagctctcttacagaagcagggtatccattccccagtgacttggagctcacagaagaggagctgagggaaccagacgacaactccctccagtcctctgaaatgaaatgttacagagccagaaactatggca
aatcaacaaccaagagaagtcctggcgaagggccccctcctcaatgttagccaccttattgcccctagcattaagaacaacccccacgtggcctcagtggtcaggagatctacaatttcacccgacagcatttccctttttctggacagctcggatggctggaagagcaccattcattacaac
ctctgcttctggacagcttggagaaggtgccagacagccttaaggatgaagataatgcaagacctcgtcttgcctttggaagctcactaaggaggggcaccgccgcttttggaggagactggtgtcttagcctttgctcaaggagagaatccaagatgcatgagtcagccagatt
gttgacctctcttctgactttga