

PDZD11-PDZ domain containing 11 Gene

Reference: PTXBC012996

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AIPP1; PDZK11; PDZ domain-containing protein 11; ATPase-interacting PDZ protein; PMCA-interacting single-PDZ protein; plasma membrane calcium ATPase-interacting single-PDZ protein; PDZ domain containing 11		

For research use only.

GENE DESCRIPTION	PDZ domain containing 11
GENE ID	51248
ORIGIN SPECIES	Human
ACCESSIONS	BC012996
VECTOR	pDONR223
NCBI SYMBOL	PDZD11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagccggattcctatgatgactaccgggtgtttcttgcctgctatgagaatctccagcatggattcctctcatgagaggggtacaccacccggactacaacaatgagttgacccagttctgccccgaaccatcacactgaagaagcctcctggagctcagttgggatttaacat
ccgaggaggaaaggcctccagctaggcatcttcatctcaaggtgattcctgactctgacacatagagcaggactgcaggaaggggaccaagtctagctgtgaatgatgtggatttcaagatattgagcacagcaaggctgttgagatcctgaagacagctcgtgaaatcagcatg
cgtgtgcgtcttccctacaattatcatcgccaaaaagagaggactgtgcactag