



Reference: PTXBC013000

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLI-Kruppel family member GLI1; zinc finger protein GLI1; GLI; glioma-associated oncogene 1; glioma-associated oncogene homolog 1 (zinc finger protein); oncogene GLI; GLI family zinc finger 1		

For research use only.

GENE DESCRIPTION	glioma-associated oncogene homolog 1 (zinc finger protein)
GENE ID	2735
ORIGIN SPECIES	Human
ACCESSIONS	BC013000
VECTOR	pDONR223
NCBI SYMBOL	GLI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr