

ANP32B-acidic (leucine-rich) nuclear phosphoprotein 32 family, member B Gene

Reference: PTXBC013003

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APRIL; PHAPI2; SSP29; acidic leucine-rich nuclear phosphoprotein 32 family member B; acidic (leucine-rich) nuclear phosphoprotein 32 family, member B; acidic protein rich in leucines; silver-stainable protein SSP29; acidic nuclear phosphoprotein 32 family		

For research use only.

GENE DESCRIPTION	acidic (leucine-rich) nuclear phosphoprotein 32 family, member B
GENE ID	10541
ORIGIN SPECIES	Human
ACCESSIONS	BC013003
VECTOR	pDONR223
NCBI SYMBOL	ANP32B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacatgaagaggagatccacctggagctgaggaacccgaccccgagctgttcgagaactgtcttgacaattgcaaatcaaatgatggaaaaattgagggttaacagctgaatttgaacttagagttcctcagtttaataatgtaggcttgatctcagtttcaaatctcccaa
gctgcctaaattgaaaaagcttgaactcagtgaaaatagaatcttggaggtctggacatgttagctgaaaaacttccaaatctcacacatctaaacttaagtggaaataaactgaaagatcagcaccttggaaccttggaaaaagttagaatgtctgaaaagcctggacctctttaactgtgag
gttaccacctgaatgactaccgagagagtgcttcaagctcctgccccagcttacctacttggatggctatgaccgagaggaccaggagcacctgactcagatgccgaggtggatgtgtggtgaggaaggaggagacgaagaaggagaagatgagggaagcaggacgatgag
gatggtgaagaagaggaggttggatgaagaagatgaagatgaagatgtagaaggggatgagcagcagatgaagtcagtggaggagaagaagaatttgacttgatgaagaagatgaagatgaggatgaggatgaaggaggagaaggtgggaaagggtgaaaagggaag
agagaacacatgatgaaggagaagatgattaa