

PSMB3-proteasome (prosome, macropain) subunit, beta type, 3 Gene

Reference: PTXBC013008

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HC10-II; proteasome subunit beta type-3; proteasome (prosome, macropain) subunit, beta type, 3; proteasome chain 13; proteasome component C10-II; proteasome theta chain; proteasome subunit beta 3		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, beta type, 3
GENE ID	5691
ORIGIN SPECIES	Human
ACCESSIONS	BC013008
VECTOR	pDONR223
NCBI SYMBOL	PSMB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctattatgtctataacggaggggcgctcatggccatgaagggaagaactgtgtggccatgctgcagacaggcgcttcgggatccaggccagatgggtaccacgggacttcagaagatcttccatgggtgaccggctgtacatcggtctggccgggctgcactgacgtc
cagacagtgtccagcgccctcaagttccggctgaacctgtatgagttgaaggaaaggtcggcagatcaaaacctataccctcatgagcatgggtggccaacctctgtatgagaaacggtttggcccttactacactgagccatgaltggcgggttgaccggaagaccttaagccctcatit
gtctctagacctcatcggtgccccatgggtgactgatgacttgggtcagtgccacctgcgccgaacaaatgtacgggaatgtgtgagtcctctgggagcccaacatggatccggatcacctgttgaacacatcccaagccatgctgaatgctgtggaccgggatgcatgtcaggca
tgggagtcattgtccacatcatcgagaaggacaaaatcaccaccaggacactgaaggccgaatggactaa