

NUMA1-nuclear mitotic apparatus protein 1 Gene

Reference: PTXBC013023

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NMP-22; NUMA; nuclear mitotic apparatus protein 1; SP-H antigen; centrophilin stabilizes mitotic spindle in mitotic cells; nuclear matrix protein-22; structural nuclear protein		

For research use only.

GENE DESCRIPTION	nuclear mitotic apparatus protein 1
GENE ID	4926
ORIGIN SPECIES	Human
ACCESSIONS	BC013023
VECTOR	pDONR223
NCBI SYMBOL	NUMA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1