

RAB3C-RAB3C, member RAS oncogene family Gene

Reference: PTXBC013033

Product type: DNA Storage Condition: -20°C
Size/Quantity: 2 ug Delivery Condition: +4°C, Blue Ice
Synonyms: RAB3C, member RAS oncogene family; ras-related protein Rab-3C
For research use only.

GENE DESCRIPTION	RAB3C, member RAS oncogene family
GENE ID	115827
ORIGIN SPECIES	Human
ACCESSIONS	BC013033
VECTOR	pDONR223
NCBI SYMBOL	RAB3C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagacacgaagcggccatgcagatggcctctgcccaagatgccaggtacgccagaaagactcctctgacagaaacttgaactacatgccaattactcaccatggcaatagcagtggtgggaaaacatctttctattcgttatgcagatgactcctttacatctgattcgtcagcaca
gtgggatcgattcaagtaaaaactgtattcaaaaatgaaaagagaatcaagcttcagatttgggacacagcaggccaggaagatacaggactaccaccagcctattatcgtggagccatgggctttatttaatgatgacattacaaatgaagaatcctcaatgcagtacaagattgg
tcaactcaaatcaaacatactcttgggacaatgcccaagttattctgttgggaacaagtgtgacatggaagacgagcgggtcactcaactgagcaggtcaacatttaggagaacagctgggtttgagttttgaaacaagtgcccaaggacaacattaatgcaagcagacattgagc
gccttgggatcactcgcgacaaaatgacagagttggagactgatcctgccatcactgctgcaagcagaacacgagactcaaggaaactcctcctccaccgagcccaactgtgcctgctag