



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

REEP4-receptor accessory protein 4 Gene

Reference: PTXBC013048

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

C8orf20; PP432; Yip2c; receptor expression-enhancing protein 4; receptor accessory protein 4

For research use only.

GENE DESCRIPTION	receptor accessory protein 4
GENE ID	80346
ORIGIN SPECIES	Human
ACCESSIONS	BC013048
VECTOR	pDONR223
NCBI SYMBOL	REEP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]