

RRAS2-related RAS viral (r-ras) oncogene homolog 2 Gene

Reference: PTXBC013106

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TC21; ras-related protein R-Ras2; ras-like protein TC21; teratocarcinoma oncogene; related RAS viral (r-ras) oncogene homolog 2		

For research use only.

GENE DESCRIPTION	related RAS viral (r-ras) oncogene homolog 2
GENE ID	22800
ORIGIN SPECIES	Human
ACCESSIONS	BC013106
VECTOR	pDONR223
NCBI SYMBOL	RRAS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcccggcgtggcgggacgctccggcaggagaagtaccggtcgtgggtggtcggcggggcggcgtggcgaagtcggcgtcaccatccagttccagtcctatttgaacggattatgacaaacattgaagattctacaaaagcagtggtgatagatgacagagcagcccggctagatatttggatacagcaggacaagaagatttggagccatgagagaacagtatatgaggactggcgaaggcttcctgttggcttttcagtcacagatagaggcagtttgaagaaatctataagttcaaaagacagattctcagagtaaggatcgtgatgagttccaatgattttaattgtaataaagcagatctggatcatcaaaagacaggttaacacaggaagaaggacaacagtttagcacggcagcttaaggtaacatacatggaggcatcagcaaaagattaggatgaatgtagatcaagcttccatgaactgtccgggttatcaggaaattcaagagcaggaatgtcctctcaccagaaccaacacggaaagaaaaagacaagaaggctgccattgtgtcatttttag