

WDR42A-WD repeat domain 42A Gene

Reference: PTXBC013107

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WDR42A; GAN2; H326; DDB1- and CUL4-associated factor 8; WD repeat domain 42A; WD repeat-containing protein 42A; DDB1 and CUL4 associated factor 8		

For research use only.

GENE DESCRIPTION	WD repeat domain 42A
GENE ID	50717
ORIGIN SPECIES	Human
ACCESSIONS	BC013107
VECTOR	pDONR223
NCBI SYMBOL	WDR42A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccagcaaaaggagcagcacagatggcagaacagacttagctaatggaagcctgtctagcagtcagagagatgtctggagctgaagagggaggagagacatcctcagccattgaagtggagcctcagacctgagttgagcttgactggggatgatggtggcccaaccgc
accagcacagaaagtcgagggcacagagacagagagctcaggtgaagataaggactctgacagcatggaggacactggtcattactccattaatgatgaaatcaggtccatgacctcagaggaagaggagagggaagaagaggaggaagaagagcagcctcgccgct
gtacagcgcaagcgggctaaccgtgaccaggactcagatgatgagcgggcccagagagactgggtgtctcagaacatcagctctaccccgacctcgtggaagccctccctgcccttcgggagcgggagctgggtcaagtgcctgtgtctatgagcctgtgggcaa
gagtcctttgtcagcgtttccgctgagcagtggtgagggccatactggtgtgtgtaataaccctgacatttaaccagcgcggcacctggctggccagtgccagcgaatgacctgaaggtggtgtgtgggtgggtacggcgccagccagctgagccttgagagtgcccaaaag
taattgtgtccaggtgaggcaaggagcagataatagcaactgagagaatcaggcatgagttagagaaatcgaactgcagcatggattgggagctgattgaggtattgttag