

IL1A-interleukin 1, alpha Gene

Reference: PTXBC013142

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-1A; IL1; IL1-ALPHA; IL1F1; interleukin-1 alpha; IL-1 alpha; hematopoietin-1; preinterleukin 1 alpha; pro-interleukin-1-alpha; interleukin 1 alpha		

For research use only.

GENE DESCRIPTION	interleukin 1, alpha
GENE ID	3552
ORIGIN SPECIES	Human
ACCESSIONS	BC013142
VECTOR	pDONR223
NCBI SYMBOL	IL1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaagttccagacatgtttgaagacctgaagaactgttacagtgaagaatgaagaagacagttcctcattgatcatctgtctctgaatcagaatcctctatcatgtaagctatggccactccatgaaggctcatggatcaatctgtgtctctgagtatcttgaacctctaaacatc
caagcttacctcaaggagagcatggtggttagtagcaaccaacgggaagggtctgaagaagagacggttgagtttaagccaatccatcactgatgatgacctggaggccatcgccaatgactcagaggaagaatcatcaagcctagggtcagcaccttttagcttcctgagcaatgtgaaat
acaactttatgaggatcatcaaatacgaattacatcctgaatgacgccctcaatcaagataattcagccaatgatcagttacctcagggctgtgcttaccataatctggatgaagcagtgaaatttgacatgggtgcttataagtcataaaggatgctaaaaattaccgtgattctaaagaatc
tcaaaaactcaattgtatgtgactgcccaagatgaagaccaaccagtgctgctgaaggagatgcctgagataccaaaaccatcacaggtagtgagaccaacctctctcttctgggaaactcacggcactaagaactatttcacatcagttgcccatcaaacctgtttattgcccaagca
agactactgggtgtgcttggcaggggggccacccctatcactgactttcagatactggaaaaccaggcgtag