

MRPL32-mitochondrial ribosomal protein L32 Gene

Reference: PTXBC013147

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC283; L32mt; MRP-L32; bMRP-59b; 39S ribosomal protein L32, mitochondrial; mitochondrial ribosomal protein L32		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L32
GENE ID	64983
ORIGIN SPECIES	Human
ACCESSIONS	BC013147
VECTOR	pDONR223
NCBI SYMBOL	MRPL32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtggccatgctgtcttgggtgttgcgggtgctgctggccgggagtgcttcgaaactactgggagcactgctacggaagcttcgcagagccggcgggcttccagtcctcgtggggaccagcattagcagtagggccagccatgtttacagagccagcaaat
gataccagtggaagtaaaagagaattccagccttttggacagfatctttggatggcagctcccaaaatagacgcacattgaagttaaccgggtgtaggagaagaatccgcagaagcttataaagttaagaacaacatagacgtttgtcctgaatgtgttcacctgaacagaaacatgtcc
tttgtgcctactgctatgaaaagggtgtgcaaggagactgcagaatcagacgacagatagggaagcaagaaggggcccttttaaggctccaccatagagactgtgtgtgtgtacacaggagagacaccgtctgaacaagatcagggaaggatcattgaacgagacagaagcg
accatcctggttcaccagaattga