

CRYZL1-crystallin, zeta (quinone reductase)-like 1 Gene

Reference: PTXBC013155

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4P11; QOH-1; quinone oxidoreductase-like protein 1; crystallin, zeta (quinone reductase)-like 1; protein 4P11; quinone oxidoreductase homolog 1; quinone reductase-like 1; zeta-crystallin homolog; crystallin zeta like 1		

For research use only.

GENE DESCRIPTION	crystallin, zeta (quinone reductase)-like 1
GENE ID	9946
ORIGIN SPECIES	Human
ACCESSIONS	BC013155
VECTOR	pDONR223
NCBI SYMBOL	CRYZL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaggcttattccaacagagttccacagatgaagaaataacattgtattcaagaaaaggaagatcttctgttacagaggataacttgtgaaactcaagttaaagcttgctctgagccagataatacaaaagcttctggcagaaatgaagatgaaaaggatttattctgttggg
agagaaattgctggaattgtattagattgtggaagcaaggatcattcttcaaccagatgatgaagtagtgaagtacttcccggctgggtgctgtggtcagctctgtaatccagcacgttgggagccaaagtggtggcgaatcacctga