

SELK-selenoprotein K Gene

Reference: PTXBC013162

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SELK; HSPC030; HSPC297		
For research use only.			

GENE DESCRIPTION	selenoprotein K
GENE ID	58515
ORIGIN SPECIES	Human
ACCESSIONS	BC013162
VECTOR	pDONR223
NCBI SYMBOL	SELK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtttacatctcgaacggacaagtgttgacagccggagtcagtcctcatggagattatctttgataacagattctctggggaatagctgagttgtgtttttttcaaaactctgttcagcaagatgtgaaaaaagaagaagctatgaaactcatctgattccagatatgatgatga
agagggccaccaggaaacctccccgaagaatgggtagaatcaatcatctgcgtggccctagtcctccctccaatggctggtgagtaggaaggttaa