

DCTN6-dynactin 6 Gene

Reference: PTXBC013175

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WS-3; WS3; p27; dynactin subunit 6; dynactin 6; dynactin subunit p27; novel RGD-containing protein; protein WS-3		

For research use only.

GENE DESCRIPTION	dynactin 6
GENE ID	10671
ORIGIN SPECIES	Human
ACCESSIONS	BC013175
VECTOR	pDONR223
NCBI SYMBOL	DCTN6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagaagactcaaaagagtgtgaagattgctcctggagcagttgtatgtgtagaaagtgaatcagaggagatgtaactatcgacacagcagtgatccacctaagcaagaattatgcggagccgggccaatagtgattggcgaagggaacctaataagaacagggccttatcataaatgcttaccagataataatcactcctgacactgaagatccagaacaaaacatgatcattggcaccaataatgtgttgaagttggctgtattcccaagccatgaagatgggagataataatgcatgaaatcaaaagcatatgtaggcagaaatgtaattgacaagtggtgcacattgggctgttgcacctaataacattgaagtcacccctgagaatacgggtgatctatgggtgcagactgcctctgctgggtgcagactgagcgaccgcagcccagacactacagctggattctctgatgaaaaatcttgccaaattaccaccacctaagaagactatgaaaggaagctcaactccagtaaaactaa