



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SFRS16-splicing factor, arginine/serine-rich 16 Gene

Reference: PTXBC013178

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

SFRS16; CLASP; SWAP2; CLK4-associating serine/arginine rich protein; Clk4 associating SR-related protein; splicing factor, arginine/serine-rich 16 (suppressor-of-white-apricot homolog, *Drosophila*); suppressor of white apricot homolog 2; CLK4 associating s

For research use only.

GENE DESCRIPTION	splicing factor, arginine/serine-rich 16
GENE ID	11129
ORIGIN SPECIES	Human
ACCESSIONS	BC013178
VECTOR	pDONR223
NCBI SYMBOL	SFRS16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr