

TRIM65-tripartite motif-containing 65 Gene

Reference: PTXBC013181

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4732463G12Rik; tripartite motif-containing protein 65; tripartite motif containing 65		
For research use only.			

GENE DESCRIPTION	tripartite motif-containing 65
GENE ID	201292
ORIGIN SPECIES	Human
ACCESSIONS	BC013181
VECTOR	pDONR223
NCBI SYMBOL	TRIM65
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccagagcttcaggccggcaccactactgggaggtgcgcgcgtcagaccactcggtagactggcgctctctacccgcaactgccacggcgcaggtggggcccccacacagacaacattggccggggaccctgctcctggggcctctgcgtccaggaggacagcctcc
aggcctggcacaacggggaagcccagcgcctcccaggggtgcaggcggcctcctggcgtggttggacctggcctcaggtgcctcactctctacagcctggagccccagaccagccctgtacacctccatgcctcttcaaccagccctcaccctcgtctctgctcctcg
agggtag