

SNCA-synuclein, alpha (non A4 component of amyloid precursor) Gene

Reference: PTXBC013293

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NACP; PARK1; PARK4; PD1; alpha-synuclein; non A-beta component of AD amyloid; synuclein alpha-140; synuclein, alpha (non A4 component of amyloid precursor); synuclein alpha		

For research use only.

GENE DESCRIPTION	synuclein, alpha (non A4 component of amyloid precursor)
GENE ID	6622
ORIGIN SPECIES	Human
ACCESSIONS	BC013293
VECTOR	pDONR223
NCBI SYMBOL	SNCA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgtattcatgaaaggactttcaaggccaaggaggaggtgtgtgctgctgagaaaaacaaacaggggtgtgcagaagcagcaggaagacaaaagggtgttctctatgtaggtctcaaaacaaaggaggaggtgtgtgcacagtggtgagaagacca
aagagcaagtgcacaaatgttgaggagcagtggtgacgggtgtgacagcagtagccagaagacagtgaggaggagcaggagcattgcagcagccactggtttgtcaaaaaggaccagttggcagaagaatgaagaaggagcccccacaggaaggaaattctggaagatagcctgt
ggatcctgacaaatgaggcttatgaaatgcctctgaggaagggtatcaagactacgaacctgaagcctaa